

USER GUIDE

Ammonia Water Test Kit






Be sure to remember

1. Make sure the vial is clean or wash the test vial properly before the test.
2. To get maximum accuracy always place the test sample vial against white background.
3. Test the samples immediately after collecting it from the source.

Directions for Testing Ammonia Levels

1. Take 4ML water in vial
2. Add 5 drops **Ammonia Reagent 1**
3. Add 5 drops **Ammonia Reagent 2**
4. Add 5 drops **Ammonia Reagent 3**
5. Caped the Vial and shake vigorously for 5 seconds
6. Wait for few (10-15) minutes
7. Match the vial colour only from provided

TRM Ammonia (Amonique) Colour Card



Wait for few (10-15) minutes

Step - 6



Note: DO NOT hold finger over the open end of the tube, as this may affect the test results. Wait 5 minutes for the colour to develop. See the test results by matching the test solution against the **Ammonia Test Colour Chart**.